



MAP3K1 Antibody

Product Code	CSB-PA442614
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13233
Immunogen	Synthesized peptide derived from MEK Kinase-1 . at AA range: 1340-1420
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	mitogen-activated protein kinase kinase kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.