




PIK3C3 Antibody

Product Code	CSB-PA381950
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEB9
Immunogen	Synthetic Peptide of PI 3 Kinase Class 3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Cow
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	phosphoinositide-3-kinase, class 3
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3C3

Image



Western Blot analysis of 293T 3T3 cells using PI 3 Kinase Class 3 Polyclonal Antibody diluted at 1:1500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.