




TRAF1 Antibody

Product Code	CSB-PA118746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13077
Immunogen	Synthesized peptide derived from TRAF1 at AA range: 191-240
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	TNF receptor-associated factor 1 (Epstein-Barr virus-induced protein 6)
Immunogen Species	Homo sapiens (Human)
Target Names	TRAF1

Image



Western Blot analysis of HEPG2 cells using TRAF1 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.