




Pan myristoylation Antibody

Product Code	CSB-PA885152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Synthetic myristoylated peptide
Raised In	Rabbit
Species Reactivity	ALL
Tested Applications	WB,ELISA;WB:1:500-10000,ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Image



Western Blot analysis of KB cells using Antibody diluted at 500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.