





HSPA8 Antibody

Product Code	CSB-PA575152
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11142
Immunogen	Recombinant Protein of Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Immunogen Species	Homo sapiens (Human)
Target Names	HSPA8


Image



Immunohistochemical analysis of paraffin-embedded Human Testis Tissue using HSC70 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Testis Tissue using HSC70 Rabbit pAb diluted at 1:200



Western blot analysis of 1) HeLa Cell Lysate, 2) Mouse Brain Tissue Lysate, 3) Rat Brain Tissue Lysate using HSC70 Rabbit pAb diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.