





HSP90AA1 Antibody

Product Code	CSB-PA330916
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07900
Immunogen	Recombinant Protein of HSP90 α
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC,ELISA;IHC:1:50-200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Heat shock protein HSP 90-alpha (Heat shock 86 kDa) (HSP 86) (HSP86) (Renal carcinoma antigen NY-REN-38)
Immunogen Species	Homo sapiens (Human)
Target Names	HSP90AA1

Image



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using HSP90 α Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using HSP90 α Rabbit pAb diluted at 1:200