





MAP2K2 Antibody

Product Code	CSB-PA101222
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36507
Immunogen	Recombinant Protein of MEK2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Dual specificity mitogen-activated protein kinase kinase 2 (MAP kinase kinase 2) (MAPKK 2) (EC 2.7.12.2) (ERK activator kinase 2) (MAPK/ERK kinase 2) (MEK 2)
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K2

Image



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using MEK2 Rabbit pAb diluted at 1:200



Western blot analysis of 1) HeLa Cell Lysate, 2) 293t Cell Lysate using MEK2 Rabbit pAb diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.