





# PTK2 Antibody

<b>Product Code</b>	CSB-PA147655
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q05397
<b>Immunogen</b>	Recombinant Protein of FAK
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	IHC,ELISA;IHC:1:50-200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (Protein phosphatase 1 regulatory subunit 71) (PPP1R71) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PTK2

## Image



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using FAK Rabbit pAb diluted at 1:500



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using FAK Rabbit pAb diluted at 1:500

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.