




Phospho-MDM2 (Tyr394) Antibody

Product Code	CSB-PA455127
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00987
Immunogen	Synthetic peptide from Human protein at AA range: 330-410
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB: 1:500-2000, ELISA: 1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	E3 ubiquitin-protein ligase Mdm2 (EC 6.3.2.-) (Double minute 2 protein) (Hdm2) (Oncoprotein Mdm2) (p53-binding protein Mdm2)
Immunogen Species	Homo sapiens (Human)
Target Names	MDM2

Image



Western blot analysis of HELA-UV Cell Lysate using antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.