




Phospho-PRC1 (Thr481) Antibody

Product Code	CSB-PA399087
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43663
Immunogen	Synthetic Phospho peptide from Human protein at AA range: 481
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Protein regulator of cytokinesis 1
Immunogen Species	Homo sapiens (Human)
Target Names	PRC1

Image



Western blot analysis of 293T VEC lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.