




PF4 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA958757 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02776 |
| Immunogen | Synthetic peptide from Human protein at AA range: 51-100 |
| Raised In | Rabbit |
| Species Reactivity | Human,Rat |
| Tested Applications | IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Alias | Platelet factor 4 (PF-4) (C-X-C motif chemokine 4) (Iroplact) (Oncostatin-A) [Cleaved into: Platelet factor 4, short form] |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PF4 |

Image



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200