





IFNAR1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA107635 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17181 |
| Immunogen | Synthetic peptide from Human protein at AA range: 300-350 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Alias | Interferon alpha/beta receptor 1 (IFN-R-1) (IFN-alpha/beta receptor 1) (Cytokine receptor class-II member 1) (Cytokine receptor family 2 member 1) (CRF2-1) (Type I interferon receptor 1) |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IFNAR1 |

Image



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.