




LGALS3BP Antibody

Product Code	CSB-PA667994
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08380
Immunogen	Synthetic peptide from Human protein at AA range: 100-160
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LGALS3BP

Image



Western blot analysis of human blood lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.