




SEMA7A Antibody

Product Code	CSB-PA613575
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75326
Immunogen	Synthetic peptide from Human protein at AA range: 101-150
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Semaphorin-7A (CDw108) (JMH blood group antigen) (John-Milton-Hargen human blood group Ag) (Semaphorin-K1) (Sema K1) (Semaphorin-L) (Sema L) (CD antigen CD108)
Immunogen Species	Homo sapiens (Human)
Target Names	SEMA7A

Image



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.