




# ZAP70 Antibody

<b>Product Code</b>	CSB-PA823016
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P43403
<b>Immunogen</b>	Synthetic peptide from Human protein at AA range: 391-440
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Tyrosine-protein kinase ZAP-70 (EC 2.7.10.2) (70 kDa zeta-chain associated protein) (Syk-related tyrosine kinase)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ZAP70

## Image



Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:200