





UBB/UBC/RPS27A/UBA52 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA374539 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0CG47/P0CG48/P62979/P62987 |
| Immunogen | Synthetic peptide from Human protein at AA range: 1-50 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | UBB/UBC/RPS27A/UBA52 |

Image



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.