




# NFATC2 Antibody

<b>Product Code</b>	CSB-PA418194
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q13469
<b>Immunogen</b>	Synthetic peptide from Human protein at AA range: 640-700
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Nuclear factor of activated T-cells, cytoplasmic 2 (NF-ATc2) (NFATc2) (NFAT pre-existing subunit) (NF-ATp) (T-cell transcription factor NFAT1)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NFATC2

## Image



Immunohistochemical analysis of paraffin-embedded Human-colon, antibody was diluted at 1:100

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.