




NAIP Antibody

Product Code	CSB-PA972562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13075
Immunogen	Synthetic peptide from Human protein at AA range: 1191-1240
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	IHC,ELISA;IHC:1:50-300,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Baculoviral IAP repeat-containing protein 1 (Neuronal apoptosis inhibitory protein)
Immunogen Species	Homo sapiens (Human)
Target Names	NAIP

Image



Immunohistochemical analysis of paraffin-embedded Human-placenta, antibody was diluted at 1:100