




NAIP Antibody

Product Code	CSB-PA972562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13075
Immunogen	Synthetic peptide from Human protein at AA range: 1191-1240
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC: 1:50-300, ELISA: 1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Baculoviral IAP repeat-containing protein 1 (Neuronal apoptosis inhibitory protein)
Immunogen Species	Homo sapiens (Human)
Target Names	NAIP

Image



Immunohistochemical analysis of paraffin-embedded Human-placenta, antibody was diluted at 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.