




ACO2 Antibody

Product Code	CSB-PA920382
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99798
Immunogen	Synthesized peptide derived from Human ACO2. at AA range: 421-470
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Alias	Aconitate hydratase, mitochondrial (Aconitase) (EC 4.2.1.3) (Citrate hydro-lyase)
Immunogen Species	Homo sapiens (Human)
Target Names	ACO2

Image



Western Blot analysis of 1,mouse-kidney
 2,mouse-heart 3,3T3 4,Hela cells using primary
 antibody diluted at 1:500(4°C overnight).
 Secondary antibody:Goat Anti-rabbit IgG IRDye
 800(diluted at 1:5000, 25°C, 1 hour)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.