




# SMARCB1 Antibody

<b>Product Code</b>	CSB-PA892944
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q12824
<b>Immunogen</b>	Synthetic peptide from AA range: 331-380
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB,ELISA;WB:1:500-10000ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	SMARCB1 BAF47 INI1 SNF5L1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SMARCB1

## Image



Western Blot analysis of HEPG2 cells using Antibody diluted at 800

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.