




KRT14/KRT17 Antibody

Product Code	CSB-PA563215
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02533/Q04695
Immunogen	Synthetic peptide from AA range: 261-310
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-10000ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	KRT14
Immunogen Species	Homo sapiens (Human)
Target Names	KRT14/KRT17

Image



Western Blot analysis of HEPG2 cells using Antibody diluted at 1000