




FBXO7 Antibody

Product Code	CSB-PA741109
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3I1
Immunogen	Synthetic peptide from AA range: 371-420
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB,ELISA;WB:1:500-10000ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	FBXO7 FBX7
Immunogen Species	Homo sapiens (Human)
Target Names	FBXO7

Image



Western Blot analysis of mouse-kidney cells using Antibody diluted at 1000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.