




SUPT20H Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA394854 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NEM7 |
| Immunogen | Synthetic peptide from AA range: 231-280 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB,ELISA;WB:1:500-10000ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Alias | FAM48A C13orf19 FP757 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | SUPT20H |

Image



Western Blot analysis of mouse-kidney cells using Antibody diluted at 1000