




CHID1 Antibody

Product Code	CSB-PA740982
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWS9
Immunogen	Synthetic peptide from AA range: 81-130
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-10000ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	CHID1 GL008 PSEC0104 SB139
Immunogen Species	Homo sapiens (Human)
Target Names	CHID1

Image



Western Blot analysis of 3T3 cells using
Antibody diluted at 500