




CHID1 Antibody

Product Code	CSB-PA740982
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWS9
Immunogen	Synthetic peptide from AA range: 81-130
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB,ELISA;WB:1:500-10000ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	CHID1 GL008 PSEC0104 SB139
Immunogen Species	Homo sapiens (Human)
Target Names	CHID1

Image



Western Blot analysis of 3T3 cells using Antibody diluted at 500

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.