




# TRPV4 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA160091  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q9HBA0  |
| <b>Immunogen</b>           | Synthetic peptide from AA range: 461-510  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human,Mouse,Rat   |
| <b>Tested Applications</b> | WB,ELISA;WB:1:500-10000ELISA:1:10000  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                   |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| <b>Alias</b>               | TRPV4 VRL2 VROAC  |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | TRPV4   |

## Image



Western Blot analysis of HEPG2 cells using Antibody diluted at 500

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.