





MAP2K3 Antibody

Product Code	CSB-PA211696
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46734
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Dual specificity mitogen-activated protein kinase kinase 3 (MAP kinase kinase 3) (MAPKK 3) (EC 2.7.12.2) (MAPK/ERK kinase 3) (MEK 3) (Stress-activated protein kinase kinase 2) (SAPK kinase 2) (SAPKK-2) (SAPKK2)
Immunogen Species	Homo sapiens (Human)
Target Names	MAP2K3

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using MEK3 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using MEK3 Rabbit pAb diluted at 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.