





ATG5 Antibody

Product Code	CSB-PA334611
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1Y0
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)
Immunogen Species	Homo sapiens (Human)
Target Names	ATG5

Image



Immunohistochemical analysis of paraffin-embedded Human Hepatocarcinoma Tissue using ATG5 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Ovarian Carcinoma Tissue using ATG5 Rabbit pAb diluted at 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.