





ATG7 Antibody

Product Code	CSB-PA443647
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95352
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Ubiquitin-like modifier-activating enzyme ATG7 (ATG12-activating enzyme E1 ATG7) (Autophagy-related protein 7) (APG7-like) (hAGP7) (Ubiquitin-activating enzyme E1-like protein)
Immunogen Species	Homo sapiens (Human)
Target Names	ATG7

Image



Immunohistochemical analysis of paraffin-embedded Human Hepatocarcinoma Tissue using ATG7 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ATG7 Rabbit pAb diluted at 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.