





# ATG13 Antibody

<b>Product Code</b>	CSB-PA560481
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O75143
<b>Immunogen</b>	Recombinant Protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Autophagy-related protein 13
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ATG13

## Image



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using ATG13 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using ATG13 Rabbit pAb diluted at 1:200

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.