





ATG13 Antibody

Product Code	CSB-PA560481
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75143
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Autophagy-related protein 13
Immunogen Species	Homo sapiens (Human)
Target Names	ATG13

Image



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using ATG13 Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using ATG13 Rabbit pAb diluted at 1:200