





ABL1 Antibody

Product Code	CSB-PA959474
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00519
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Tyrosine-protein kinase ABL1 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 1) (Abelson tyrosine-protein kinase 1) (Proto-oncogene c-Abl) (p150)
Immunogen Species	Homo sapiens (Human)
Target Names	ABL1

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using c-Abl Rabbit pAb diluted at 1:200



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using c-Abl Rabbit pAb diluted at 1:200

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.