





NBR1 Antibody

Product Code	CSB-PA286085
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14596
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	IHC, ELISA; IHC:1:50-1:200, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Next to BRCA1 gene 1 protein (Cell migration-inducing gene 19 protein) (Membrane component chromosome 17 surface marker 2) (Neighbor of BRCA1 gene 1 protein) (Protein 1A1-3B)
Immunogen Species	Homo sapiens (Human)
Target Names	NBR1

Image



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using NBR1 Rabbit pAb diluted at 1:500



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using NBR1 Rabbit pAb diluted at 1:500