




Histone H2A.Z Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA393326 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0C0S5 |
| Immunogen | Synthetic peptide from AA range: 1-30 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB,ELISA;WB:1:500-2000,ELISA:1:10000-20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Alias | Histone H2A.Z (H2A/z) |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | H2AFZ |

Image



Western blot analysis of HELA Cell Lysate using antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.