




COL12A1 Antibody

Product Code	CSB-PA669302
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99715
Immunogen	Synthetic peptide from Human protein at AA range: 1268-1341
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC,ELISA;IHC:1:50-1:300,ELISA:1:5000-1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Collagen alpha-1(XII) chain
Immunogen Species	Homo sapiens (Human)
Target Names	COL12A1

Image



Immunohistochemical analysis of paraffin-embedded human-muscle, antibody was diluted at 1:200