




PLK2 Antibody

Product Code	CSB-PA745425
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYY3
Immunogen	Synthetic peptide from Human protein at AA range: 1268-1341
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	IHC,ELISA;IHC:1:50-1:300,ELISA:1:5000-1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Alias	Serine/threonine-protein kinase PLK2 (EC 2.7.11.21) (Polo-like kinase 2) (PLK-2) (hPlk2) (Serine/threonine-protein kinase SNK) (hSNK) (Serum-inducible kinase)
Immunogen Species	Homo sapiens (Human)
Target Names	PLK2

Image



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.