




# PLK2 Antibody

<b>Product Code</b>	CSB-PA745425
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9NYY3
<b>Immunogen</b>	Synthetic peptide from Human protein at AA range: 1268-1341
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Rat
<b>Tested Applications</b>	IHC,ELISA;IHC:1:50-1:300,ELISA:1:5000-1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Serine/threonine-protein kinase PLK2 (EC 2.7.11.21) (Polo-like kinase 2) (PLK-2) (hPlk2) (Serine/threonine-protein kinase SNK) (hSNK) (Serum-inducible kinase)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PLK2

## Image



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200