




# Phospho-TP53 (Ser366) Antibody

<b>Product Code</b>	CSB-PA847762
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04637
<b>Immunogen</b>	Synthetic phospho-peptide from Human protein at phospho Ser366
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB,ELISA;WB:1:500-1:2000,ELISA:1:5000-1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	TP53

## Image



Western blot analysis of HELA Cell Lysate, antibody was diluted at 1:1000

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.