




# CASP3 Antibody

<b>Product Code</b>	CSB-PA140280
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P42574
<b>Immunogen</b>	Synthetic peptide from Human protein at AA range: 1268-1341
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	IHC,ELISA;IHC:1:50-1:300,ELISA:1:5000-1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300;Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Alias</b>	Caspase-3 (CASP-3) (EC 3.4.22.56) (Apopain) (Cysteine protease CPP32) (CPP-32) (Protein Yama) (SREBP cleavage activity 1) (SCA-1) [Cleaved into: Caspase-3 subunit p17; Caspase-3 subunit p12]
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CASP3

## Image



Immunohistochemical analysis of paraffin-embedded human-tonsils, antibody was diluted at 1:200