




# Rubisco (Large Chain) Monoclonal Antibody

<b>Product Code</b>	CSB-MA141458
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Research Use</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Protein
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Plants
<b>Tested Applications</b>	ELISA,WB;Recommended dilution:WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Product Type</b>	Monoclonal Antibody

## Image



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at 1)1:2,000?2)1:5,000