




Rubisco (Large Chain) Monoclonal Antibody

Product Code	CSB-MA141458
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Protein
Raised In	Mouse
Species Reactivity	Plants
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Immunogen Species	Plants

Image



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at 1)1:2,000 2)1:5,000