




Rubisco (Large Chain) Monoclonal Antibody

Product Code	CSB-MA271230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Protein
Raised In	Mouse
Species Reactivity	Plants
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody

Image



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at (1)1:2,000;(2)1:5,000