





TBP Monoclonal Antibody

Product Code	CSB-MA437236
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20226
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;Recommended dilution:IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Alias	TBP; GTF2D1; TF2D; TFIID; TATA-box-binding protein; TATA sequence-binding protein; TATA-binding factor; TATA-box factor; Transcription initiation factor TFIID TBP subunit
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	TBP

Image



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using TBP/TATA Binding Protein(mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using TBP/TATA Binding ProteinMouse mAb diluted at 1:200.