




S1-Tag Monoclonal Antibody

Product Code	CSB-MA649669
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Clone No.	15B5

Image

Western blot analysis of Recombinant S1-Tag Protein with S1-Tag Mouse mAb diluted at 1:5,000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.