



Digoxin Monoclonal Antibody

Product Code	CSB-MA291797
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Digoxin
Raised In	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Clone No.	28C10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.