



# Digoxin Monoclonal Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-MA291797   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Research Use</b>        | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.                  |
| <b>Immunogen</b>           | Digoxin  |
| <b>Raised In</b>           | Mouse  |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Purification Method</b> | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| <b>Product Type</b>        | Monoclonal Antibody  |