





BAX Monoclonal Antibody

Product Code	CSB-MA183751
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07812
Immunogen	Synthetic Peptide of Bax at AA range of 30-110
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Alias	BAX
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAX
Clone No.	13H6


Image




Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat Testis Tissue using Bax Mouse mAb diluted at 1:200.



Western blot analysis of 1) Hela Cell Lysate, 2) PC12 Cell Lysate using Bax Mouse mAb diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.