




Luciferase Monoclonal Antibody

Product Code	CSB-MA315104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08659
Immunogen	Recombinant Protein of Luciferase of Luciferase
Raised In	Mouse
Species Reactivity	Firefly
Tested Applications	ELISA,WB,IHC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Product Type	Monoclonal Antibody
Immunogen Species	Photinus pyralis (Common eastern firefly) (Lampyris pyralis)

Image



Western Blot analysis of Luciferase protein using antibody diluted at 1:1000