





CCNB1 Monoclonal Antibody

Product Code	CSB-MA820649
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Research Use	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14635
Immunogen	Synthetic Peptide of Cyclin B1 at AA range of 60-140
Raised In	Mouse
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC; Recommended dilution: IHC: 1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Alias	CCNB1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	28D1

Image



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Cyclin B1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Cyclin B1 Mouse mAb diluted at 1:200.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.